



100 W. Industry Ct. Deer Park, NY 11729
631-333-4100 www.SteriTatt.com

Ultra Violet Pink Ink Tips:

Keep in mind the UV Pink ink is a thinner ink and is a bit more challenging to use.

Key points to educate Therapists when introducing:

Shake UV cartridges well prior to use

Attach large gauge needle, 18g or 20g.

The Dropping technique is not recommended for UV inks. The Vaccination technique is more appropriate and will ensure getting the ink under the skin.

When first using the UV ink it is important to get used to its performance so we always suggest to use on "INSIGNIFICANT POINTS", superior border, a corner mark etc. Once the therapists are successful and comfortable with the viscosity and the adherence to the skin, at that point they can move on to using UV ink as the main ink for ISO and all points. This is good practice to avoid any issues.

Keep in mind when using the larger gauge needle the UV component in the ink sometimes takes a few days before you can see it shimmer. This is due to the skin trauma overshadowing the mark.

Please make sure to explain this in detail. I have found that it is very easy to get caught up in the "coolness" of the UV Pink and how nice it shines under the black light. So the sale is easy however if we do not educate properly and do not give realistic expectations it will be short lived.

I have had therapists say "it does not work", "I cannot see it". I attribute that to inexperience, poor technique or not enough education on the proper uses and challenges.

Please feel free to contact me with any questions. If end users want to discuss give me a heads up and then send them my email james@steritatt.com.